

# **HOMOCYSTEINE 3 REAGENT ENZYMATIC ASSAY**

**Three Reagent Liquid Stable**

Diazyme's Homocysteine 3 Reagent Enzymatic Assay features convenient ready to use reagent, calibrators and controls for the quantitative determination of total L-homocysteine in serum or plasma. Diazyme's proprietary Enzyme Cycling methodology is an excellent choice for cost conscious laboratories of all sizes due to a wide variety of instrument specific packaging options. The assay requires minimal patient sample and provides fast, accurate and precise results. A wide variety of reliable instrument parameters means the assay is readily available for installation on most automated clinical chemistry analyzers.

## **DIAZYME HOMOCYSTEINE 3 REAGENT ASSAY ADVANTAGES**

- Award winning Homocysteine recognized by the American Association of Clinical Chemistry (AACC) for outstanding contribution to scientific research
- Innovative enzyme cycling based technology for accurate and reliable results
- Excellent correlation to HPLC and immunochemical methods
- No "carry over" issues with iron or lipase reagents
- Test renal patients with confidence since there is no interference from cystathionine which affects some other less specific methods
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation saving time and reducing sample handling

## **REGULATORY STATUS**


510(k) Cleared

Health Canada Registered



## **AVAILABLE INSTRUMENT SPECIFIC PACKAGING**

- **Roche**  
- Modular P  
- Integra  
- Cobas  
- Hitachi
- **Beckman**  
- Synchron
- **Siemens**  
- Dimension



# HOMOCYSTEIN 3 REAGENT ENZYMATIC ASSAY


Three Reagent Liquid Stable



## ASSAY SPECIFICATIONS

|                                   |   |
|-----------------------------------|---|
| <b>Method</b>                     | Diazyme Patented Enzyme Cycling   |
| <b>Sample Type &amp; Volume</b>   | <ul style="list-style-type: none"><li>• Serum</li><li>• Plasma</li><li>- EDTA</li><li>- Li-heparin</li></ul> Sample Volume 18 $\mu$ L |
| <b>Method Correlation</b>         | N = 66<br>$y$ -intercept = 0.87<br>Slope = 0.98<br>$R^2$ = 0.976  |
| <b>Linear Range</b>               | up to 50 $\mu$ mol/L  |
| <b>LOD</b>                        | <1.5 $\mu$ mol/L  |
| <b>Calibration Levels</b>         | 5-Point Calibration   |
| <b>Reagent On-Board Stability</b> | Opened:<br>At least 100 days<br>(Analyzer Dependent)  |

### Homocysteine 3 Reagent Assay Procedure\*



#### \*Analyzer Dependent

Parameter questions for Enzymatic Homocysteine 3 Reagent Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email [support@diazyme.com](mailto:support@diazyme.com)

1. Vilaseca et al. Clin. Chem. 43: 690-692 (1997)
2. Faure-Delanef et al. Am. J. Hum. Genet. 60: 999-1001 (1997)

## ASSAY PRECISION

Precision studies were conducted according to the NCCLS EP-5 protocol. Within precisions (CV%) for three levels of Hcy controls are 2.2% for 7  $\mu$ M Hcy, 3.0% for 12  $\mu$ M Hcy and 1.8% for 29.5  $\mu$ M Hcy. Total imprecision for three levels of Hcy controls are 4.1% for 7  $\mu$ M Hcy, 5.9% for 12  $\mu$ M Hcy and 4.0% for 29.5  $\mu$ M Hcy.

| HCY Concentration                 | 7 $\mu$ M<br>N = 40 | 12 $\mu$ M<br>N = 80 | 29.5 $\mu$ M<br>N = 80 |
|-----------------------------------|---------------------|----------------------|------------------------|
| <b>Within-Run Imprecision CV%</b> | 2.2                 | 3.0                  | 1.8                    |
| <b>Total Imprecision CV%</b>      | 4.1                 | 5.9                  | 4.0                    |

## ASSAY INTERFERENCE

An interference study was performed by testing a serum sample spiked with varied concentrations of endogenous substances. The following substances normally present in the serum produced less than 10% deviation when tested at the stated concentrations:

|                             |             |
|-----------------------------|-------------|
| NH <sub>4</sub> Cl:         | 500 $\mu$ M |
| NaPi:                       | 1 mM        |
| NaF:                        | 1 mM        |
| Triglycerides:              | 2500 mg/dL  |
| Bilirubin:                  | 20 mg/dL    |
| Hemoglobin:                 | 1200 mg/dL  |
| *Glutathione:               | 0.5 mM      |
| Ascorbic Acid:              | 10 mM       |
| L-Cysteine:                 | 1 mM        |
| S-Adenosylmethionine (SAM): | 20 $\mu$ M  |
| **Adenosine:                | 100 $\mu$ M |
| **Cystathionine:            | 100 $\mu$ M |

\* Glutathione was originally tested at 1 mM level, the interference was +13.5%. When retested at 0.5 mM level, the interference was less than 10%.

\*\* The concentrations tested are about 5-10 times higher than the normal range of serum levels.

## REFERENCE RANGE

In most of the U.S. clinical laboratories, 15  $\mu$ mol/L is used as the cut-off value for normal level of Hcy for adults.<sup>1-2</sup> In Europe, 12  $\mu$ mol/L is used as the cut-off value. However, each laboratory is recommended to establish a range of normal values for the population in their region.

## DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064  
PO Box 85608, San Diego, CA 92186  
Tel: 858-455-4768 888-DIAZYME

[www.diazyme.com](http://www.diazyme.com) [sales@diazyme.com](mailto:sales@diazyme.com)



Tel. +41 44 456 33 33  
igz.ch igz@igz.ch

## DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland  
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366  
[sales@diazyme.de](mailto:sales@diazyme.de)

## SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park  
Shanghai, 201203, People's Republic of China  
Tel: 086-21-51320668 Fax: 086-21-51320663  
[www.lanyuanbio.com](http://www.lanyuanbio.com) [service@lanyuanbio.com](mailto:service@lanyuanbio.com)

