

HOMOCYSTEINE 2 REAGENT ENZYMATIC ASSAY

Dual Vial Liquid Stable

Diazyme's Homocysteine 2 Reagent Enzymatic Assay features convenient ready to use reagent, calibrators and controls for the quantitative determination of total L-homocysteine in serum or plasma. Diazyme's proprietary Enzyme Cycling methodology is an excellent choice for cost conscious laboratories of all sizes due to a wide variety of instrument specific packaging options. The assay requires minimal patient sample and provides fast, accurate and precise results. A wide variety of reliable instrument parameters means the assay is readily available for installation on most automated clinical chemistry analyzers.

DIAZYME HOMOCYSTEIN 2 REAGENT ASSAY ADVANTAGES

- Award winning Homocysteine recognized by the American Association of Clinical Chemistry (AACC) for outstanding contribution to scientific research
- Innovative enzyme cycling based technology for accurate and reliable results
- Excellent correlation to HPLC and immunochemical methods
- No "carry over" issues with iron or lipase reagents
- Test renal patients with confidence since there is no interference from cystathionine which affects some other less specific methods
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation saving time and reducing sample handling

REGULATORY STATUS

510(k) Cleared Health Canada Registered 

AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- **Roche**
- **Beckman**
- Hitachi
- AU Series



HOMOCYSTEIN 2 REAGENT ENZYMATIC ASSAY


Dual Vial Liquid Stable



ASSAY SPECIFICATIONS

| | |
|-----------------------------------|--|
| Method | Diazyme Patented Enzyme Cycling |
| Sample Type & Volume | <ul style="list-style-type: none">• Serum• Plasma- EDTA- Li-heparin Sample Volume 13 µL |
| Method Correlation | N = 40 y-intercept = 1.05 Slope = 0.94 R ² = 0.99 |
| Linear Range | Up to 50 µmol/L |
| LOD | 0.4 µmol/L |
| Calibration Levels | 5-Point Calibration |
| Reagent On-Board Stability | Opened: At least 60 days (Analyzer Dependent) |

Homocysteine 2 Reagent Assay Procedure*



*Analyzer Dependent

Two Reagent System

Parameter questions for Enzymatic Homocysteine 2 Reagent Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email support@diazyme.com

1. Vilaseca et al. *Clin. Chem.* 43: 690-692 (1997)
2. Faure-Delanef et al. *Am. J. Hum. Genet.* 60: 999-1001 (1997)

ASSAY PRECISION

Precision studies were conducted according to the NCCLS EP-5 protocol. Four HCY serum samples containing 7.0, 12.0, 15.6, and 29.0 µM HCY were tested.

| HCY Concentration | 7 µM | 12 µM | 15.6 µM | 29 µM |
|--|------|-------|---------|-------|
| Within-Run Imprecision CV% N = 20 | 4.5 | 1.87 | 3.04 | 2.4 |
| Total Imprecision CV% N = 30 | 5.87 | 4.88 | 5.51 | 2.57 |

ASSAY INTERFERENCE

An interference study was performed by testing a serum sample spiked with varied concentrations of endogenous substances. The following substances normally present in the serum produced less than 10% deviation when tested at the stated concentrations:

| | |
|----------------------|------------|
| Bilirubin: | 40 mg/dL |
| Triglycerides: | 1000 mg/dL |
| Hemoglobin: | 500 mg/dL |
| Bilirubin Conjugate: | 40 mg/dL |
| Ascorbic Acid: | 10 mM |
| Cystathionine: | 100 µM** |

**The concentrations tested are about 5-10 times higher than the normal range of serum levels.

REFERENCE RANGE

In most of the U.S. clinical laboratories, 15 µmol/L is used as the cut-off value for normal level of Hcy for adults.¹⁻² In Europe, 12 µmol/L is used as the cut-off value. However, each laboratory is recommended to establish a range of normal values for the population in their region.

DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064
PO Box 85608, San Diego, CA 92186
Tel: 858-455-4768 888-DIAZYME

www.diazyme.com sales@diazyme.com



Instruments
smart solutions & service
IGZ Instruments AG
Furibachstrasse 17
8107 Buchs ZH

Tel. +41 44 456 33 33
igz.ch igz@igz.ch

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Room 201, 1011 Halei Road, Zhangjiang Hi-tech Park
Shanghai, 201203, People's Republic of China
Tel: 086-21-51320668 Fax: 086-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

