

# MYOGLOBIN ASSAY

## Dual Vial Liquid Stable

Diazyme's Myoglobin Assay is an excellent cost effective cardiovascular test that is used to aid in the early detection of myocardial damage. The latex enhanced immunoturbidimetric methodology is highly precise with excellent correlations to existing commercial myoglobin tests and is designed to work on most modern high throughput clinical chemistry analyzers. The Myoglobin assay offers liquid stable reagent, calibrator, high and low end controls separately for added convenience.

### DIAZYME MYOGLOBIN ASSAY ADVANTAGES

- Fast test results for a rapid turnaround time
- Liquid stable reagent, calibrator and controls are offered separately for added convenience
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation saving time and reducing sample handling

### REGULATORY STATUS

510(k) Cleared



### AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- Roche
- Hitachi



# MYOGLOBIN ASSAY


Dual Vial Liquid Stable



## ASSAY SPECIFICATIONS

<b>Method</b>	Latex Enhanced Immunoturbidimetric
<b>Sample Type &amp; Volume</b>	<ul style="list-style-type: none"><li>• Serum</li><li>• Plasma</li><li>- Li-heparin</li><li>- K<sub>3</sub>EDTA</li></ul> Sample Volume 4 $\mu$ L
<b>Method Comparison</b>	Deming Regression: N = 66 y-intercept = -5.141 Slope = 0.959 R <sup>2</sup> = 0.9927  Regression Analysis: N = 66 y-intercept = -4.228 Slope: 0.9526 R <sup>2</sup> : 0.9855 Samples Ranged From: 16.9 – 615.9 ng/mL
<b>Assay Range</b>	13.2 – 615.9 ng/mL
<b>LOB</b>	4.4 ng/mL
<b>LOD</b>	7.2 ng/mL
<b>LOQ</b>	13.2 ng/mL
<b>Calibration Levels</b>	5-Point Calibration
<b>Reagent On-Board Stability</b>	Opened: Eight weeks on board analyzer

### Myoglobin Assay Procedure\*



#### \*Analyzer Dependent

Parameter questions for Myoglobin Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email [support@diazyme.com](mailto:support@diazyme.com)

## ASSAY PRECISION

The precision of the Diazyme Myoglobin Assay was evaluated according to CLSI EP5-A guideline. In the study, three levels of serum based controls containing approximately 66, 170, and 335 ng/mL of myoglobin, and three serum sample containing approximately 35, 150, and 414 ng/mL of myoglobin, respectively, were tested with 2 runs per day in duplicates over 20 working days. Results were calculated using the EP Evaluator software precision statistic template and summarized in the following tables:

	Control Level 1	Control Level 2	Control Level 3	Serum Level 1	Serum Level 2	Serum Level 3
<b>N</b>	80	80	80	80	80	80
<b>Mean</b>	65.97	175.8	337.0	37.78	148.6	414.3
<b>SD</b>	2.45	6.69	11.9	1.77	3.53	19.7
<b>CV%</b>	3.71%	3.87%	3.54%	4.69%	2.37%	4.80%

	Control Level 1	Control Level 2	Control Level 3	Serum Level 1	Serum Level 2	Serum Level 3
<b>N</b>	80	80	80	80	80	80
<b>Mean</b>	65.97	172.8	337.0	37.78	148.6	414.3
<b>SD</b>	3.37	7.37	14.9	1.97	5.32	21.8
<b>CV%</b>	5.10%	4.30%	4.40%	5.20%	3.58%	5.3%

## ASSAY INTERFERENCE

The following substances do not interfere with this assay at the levels tested (less than 10% bias):

Ampicillin:	up to 152 $\mu$ M	Hemoglobin:	up to 1000 mg/dL
Carbamazepine:	up to 0.13 mM	Bilirubin:	up to 40 mg/dL
Na2+-Cefoxitin:	up to 1549 $\mu$ M	Conjugated Bilirubin:	up to 40 mg/dL
Ibuprofen:	up to 2425 $\mu$ M	Triglycerides:	up to 1000 mg/dL
Cyclosporin:	up to 0.125 $\mu$ M	Intralipid:	up to 125 mg/dL
Levodopa:	up to 30.4 mM	Ascorbic acid:	up to 176 mg/dL
Methyldopa:	up to 71 $\mu$ M	Rheumatoid factor:	up to 100 IU/mL
Metronidazole:	up to 701 $\mu$ M	Heparin:	up to 1.5 IU/mL
Rifampicin:	up to 78.1 $\mu$ M	N-acetylcysteine:	up to 11.04 mM
Theophylline:	up to 222 $\mu$ M	Acetylsalicylic acid:	up to 2.78 mM
Phenylbutazone:	up to 650 $\mu$ M		
Valproic Acid,			
Sodium Salt:	up to 3.5 mM		

## DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064  
PO Box 85608, San Diego, CA 92186  
Tel: 858-455-4768 888-DIAZYME

[www.diazyme.com](http://www.diazyme.com) [sales@diazyme.com](mailto:sales@diazyme.com)



## DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland  
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366  
[sales@diazyme.de](mailto:sales@diazyme.de)

## SHANGHAI DIAZYME CO., LTD.

Room 201, 1011 Halei Road, Zhangjiang Hi-tech Park  
Shanghai, 201203, People's Republic of China  
Tel: 086-21-51320668 Fax: 086-21-51320663  
[www.lanyuanbio.com](http://www.lanyuanbio.com) [service@lanyuanbio.com](mailto:service@lanyuanbio.com)

