

N-ACETYL- β -D-GLUCOSAMINIDASE (NAG) ASSAY

Three Vial Liquid Stable

Diazyme's NAG Assay is a cost effective test that uses a "highly sophisticated and powerful colorimetric substrate" that is not affected by urine color.¹ Each kit is supplied with a calibrator set for added convenience and controls are available separately. The test offers a wide range of instrument parameters for facilitating and simplifying implementation in the laboratory. Increased NAG levels in urine are usually an early indication of renal disease and can serve as a valuable renal monitoring test in disorders such as nephritic syndrome, glomerulonephritis, drug abuse associated nephrotoxicity, diabetes-associated nephropathy, hypertension and urinary tract infections.

DIAZYME NAG ASSAY ADVANTAGES

- Fast test results (under 5.5 minutes) for a rapid turnaround time
- Liquid stable reagent, calibrator and controls are offered separately for added convenience
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation, saving time and reducing sample handling


REGULATORY STATUS

510(k) Exempt



AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- Roche
- Hitachi



NAG ASSAY


Three Vial Liquid Stable



ASSAY SPECIFICATIONS

Method	Colorimetric (Enzymatic cleavage of a colorimetric substrate)
Sample Type & Volume	• Urine Sample Volume 10 µL
Linear Range	Up to 200 U/L
LOQ	1.64 U/L
Calibration Levels	1-Point Calibration
Reagent On-Board Stability	Unopened: 24 months Opened: 1 month when stored at 2-8°C

NAG Assay Procedure*



*Analyzer Dependent

Parameter questions for NAG Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email support@diazyme.com

ASSAY PRECISION

In the study, two levels of NAG controls and one NAG urine sample containing 40.9 U/L, 124.0 U/L and 9.64 U/L NAG respectively were tested on a Hitachi 917 in one run with 20 in replicates.

Within-Run Precision:

	Sample 1	Sample 2	Sample 3
N	20	20	20
Mean	38.99	119.71	9.68
Std. Dev.	0.39	1.16	0.41
CV %	0.99%	0.97%	4.23%

ASSAY INTERFERENCE

The common urine interfering substances triglyceride, ascorbic acid, free bilirubin, and conjugated bilirubin showed no significant interference ($\geq 10\%$) up to the concentrations summarized below.

Triglyceride:	1000 mg/dL
Ascorbic Acid:	0.500 mg/dL
Bilirubin:	5 mg/dL
Bilirubin Conjugated:	5 mg/dL



IG Instruments
smart solutions & service

IGZ Instruments AG
Furtbachstrasse 17
8107 Buchs ZH

Tel. +41 44 456 33 33
igz.ch igz@igz.ch

DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064
PO Box 85608, San Diego, CA 92186
Tel: 858-455-4768 888-DIAZYME

www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Room 201, 1011 Halei Road, Zhangjiang Hi-tech Park
Shanghai, 201203, People's Republic of China
Tel: 086-21-51320668 Fax: 086-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

