

# HUMAN KAPPA (κ) AND LAMBDA (λ) FREE LIGHT CHAIN ASSAYS

## Dual Vial Liquid Stable

Diazyme’s Human Kappa (κ) and Lambda (λ) Free Light Chain assays provide accurate and reliable testing across the full analytical range and correlates well with current methods. Diazyme’s Free Light Chain Assays utilize a cost effective latex particle enhanced immunoturbidimetric methodology for use on a wide range of clinical chemistry systems. Both assays display excellent precision on the low end while its liquid stable format requires no reagent preparation saving time and reducing sample handling.

### DIAZYME HUMAN KAPPA (κ) AND LAMBDA (λ) FREE LIGHT CHAIN ASSAYS ADVANTAGES

- Fast test results (10 minutes) for a rapid turnaround time
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation saving time and reducing sample handling

### REGULATORY STATUS



510(k) Cleared

Health Canada Registered

EU:  

### AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- Roche
- Hitachi



# HUMAN KAPPA ( $\kappa$ ) AND LAMBDA ( $\lambda$ ) FREE LIGHT CHAIN ASSAYS


Dual Vial Liquid Stable Assays



## ASSAY SPECIFICATIONS

Assay	Kappa ( $\kappa$ )	Lambda ( $\lambda$ )
<b>Method</b>	Latex Enhanced Immunoturbidimetric Assay	Latex Enhanced Immunoturbidimetric Assay
<b>Sample Type &amp; Volume</b>	• Serum  Sample Volume 4.0 $\mu$ L	• Serum  Sample Volume 20 $\mu$ L
<b>Method Correlation</b>	Linear Regression N = 126 y-intercept = -2.536 Slope = 0.958 R <sup>2</sup> = 0.977  Samples range from 4.63 - 2975.8 mg/L	Linear Regression N = 126 y-intercept = 0.636 Slope = 1.054 R <sup>2</sup> = 0.972  Samples range from 6.19 to 3917.80 mg/L
<b>Linearity</b>	4.5 - 3000 mg/L	6.1 - 4000 mg/L
<b>LOB</b>	1.2 mg/L	1.7 mg/L
<b>LOD</b>	2.0 mg/L	2.9 mg/L
<b>LOQ</b>	4.5 mg/L	6.1 mg/L

### Diazyme Human Kappa ( $\kappa$ ) And Lambda ( $\lambda$ ) Free Light Chain Assay Procedures\*



\*Analyzer Dependent

For a list of validated parameters please contact  
Diazyme technical support at 858-455-4768  
or email [support@diazyme.com](mailto:support@diazyme.com)

## KAPPA ( $\kappa$ ) PRECISION

The precision of the Assay was evaluated according to CLSI EP5-A guideline. In the study, eight levels of serum specimens containing kappa FLC spanning AMR and two levels of serum based kappa FLC controls were tested with 2 runs per day with duplicates over 20 working days using multiple lots of the reagents on multiple analyzers. The precision data was analyzed according to three-way nested ANOVA and the results of mean (mg/L) and CV% are summarized below:

One lot of the reagent on three analyzers

ID	Mean n=240	Within-Run	Between Run	Between Day	Between Instrument	Total
S1	9.38	7.6%	4.9%	3.5%	5.1%	11.0%
S2	35.30	2.1%	1.5%	1.3%	2.0%	3.5%
S3	122.37	2.4%	0.9%	1.8%	0.8%	3.2%
S4	5.95	7.5%	3.0%	8.5%	N/A	11.7%
S5	15.76	2.9%	1.6%	2.1%	0.2%	3.9%
S6	25.92	1.7%	1.3%	0.5%	0.4%	2.2%
S7	139.23	1.5%	1.5%	0.7%	1.1%	2.5%
S8	2588.59	1.4%	0.6%	2.6%	1.5%	3.3%

Three lots of the reagent on one analyzer

ID	Mean n=240	Within-Run	Between Run	Between Day	Between Instrument	Total
S1	9.49	7.9%	5.7%	5.7%	2.6%	11.6%
S2	35.14	2.9%	0.8%	1.5%	0.8%	3.5%
S3	121.53	2.6%	1.5%	2.3%	N/A	3.8%
S4	5.81	6.3%	2.6%	2.7%	0.7%	7.4%
S5	15.40	2.4%	0.9%	0.6%	2.2%	3.5%
S6	25.66	1.1%	0.3%	0.6%	1.0%	1.7%
S7	138.28	1.1%	N/A	1.0%	1.0%	1.8%
S8	2588.45	1.3%	1.0%	2.5%	1.5%	3.3%

## LAMBDA ( $\lambda$ ) PRECISION

One lot of the reagent on three analyzers

ID	Mean n=240	Within-Run	Between Run	Between Day	Between Instrument	Total
S1	11.89	6.3%	4.5%	6.1%	N/A	9.9%
S2	48.65	2.1%	0.7%	1.3%	0.8%	2.7%
S3	144.55	1.6%	1.0%	1.3%	0.6%	2.4%
S4	8.31	6.5%	0.4%	3.1%	0.8%	7.3%
S5	22.27	3.0%	1.4%	2.6%	1.1%	4.3%
S6	35.09	2.1%	1.2%	1.4%	1.0%	2.9%
S7	181.13	0.8%	0.7%	1.3%	1.1%	2.0%
S8	3472.59	3.6%	N/A	2.1%	0.6%	4.2%

Three lots of the reagent on one analyzer

ID	Mean n=240	Within-Run	Between Run	Between Day	Between Instrument	Total
S1	11.78	5.9%	4.0%	7.4%	4.9%	11.4%
S2	48.95	1.3%	N/A	1.7%	0.3%	2.2%
S3	142.63	1.1%	0.8%	0.5%	0.8%	1.6%
S4	8.16	6.3%	2.3%	N/A	3.1%	7.4%
S5	21.90	2.9%	1.3%	N/A	0.6%	3.2%
S6	34.79	2.0%	0.4%	0.3%	N/A	2.1%
S7	183.02	0.6%	0.3%	0.8%	0.3%	1.1%
S8	3442.60	3.2%	N/A	2.0%	N/A	3.8%

## DIAZYME LABORATORIES, INC.

12889 Gregg Court, Poway, CA 92064  
PO Box 85608, San Diego, CA 92186  
Tel: 858-455-4768 888-DIAZYME

[www.diazyme.com](http://www.diazyme.com) [sales@diazyme.com](mailto:sales@diazyme.com)



IGZ Instruments AG  
Furtbuchsstrasse 17  
8107 Buchs ZH  
Tel. +41 44 456 33 33  
igz.ch igz@igz.ch

## DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland  
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366  
[sales@diazyme.de](mailto:sales@diazyme.de)

## SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park  
Shanghai, 201203, People's Republic of China  
Tel: 086-21-51320668 Fax: 086-21-51320663  
[www.lanyuanbio.com](http://www.lanyuanbio.com) [service@lanyuanbio.com](mailto:service@lanyuanbio.com)

