

# **25-OH VITAMIN D ASSAY**

**For Clinical Chemistry Analyzers**

**Three Vial Liquid Stable**

Diazyme is the first to offer a liquid stable and cost effective 25-OH Vitamin D Assay designed to work on fully automated clinical chemistry analyzers. The highly accurate 25-OH Vitamin D assay detects total Vitamin D (D2 + D3) and is directly traceable to NIST SRM-972 standards. The assay is highly precise and has excellent correlation to LC-MS/MS and other chemiluminescence assays. Each kit is supplied with a liquid calibrator set for added convenience. Controls are available separately.

## **DIAZYME 25-OH VITAMIN D ASSAY FOR CLINICAL CHEMISTRY ADVANTAGES**

- Time to first result (15–19 Minutes)
- Results display excellent correlations to LC-MS/MS and leading 25-OH Vitamin D Immunoassays
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation, saving time and reducing sample handling


## **REGULATORY STATUS**

510(k) Cleared      Health Canada Registered



## **AVAILABLE INSTRUMENT SPECIFIC PACKAGING**

- Roche
- Hitachi



# 25-OH VITAMIN D ASSAY

For Clinical Chemistry Analyzers


Three Vial Liquid Stable



## ASSAY SPECIFICATIONS

<b>Method</b>	FemtoQuant™ Enzyme Immunoassay
<b>Sample Type &amp; Volume</b>	<ul style="list-style-type: none"> <li>• Serum</li> <li>• Plasma</li> <li>- K<sub>3</sub>-EDTA</li> <li>- Li-heparin</li> </ul> <p>Sample Volume 20 µL</p>
<b>Method Correlation</b>	<p>N = 98</p> <p>y-intercept = -0.21 ng/mL</p> <p>Slope = 1.005</p> <p>R<sup>2</sup> = 0.968</p> <p>Sample Range: 9.5 ng/mL – 140.9 ng/mL</p>
<b>Linear Range</b>	7.6 to 147.8 ng/mL
<b>LOB</b>	2.0 ng/mL
<b>LOD</b>	3.5 ng/mL
<b>LOQ</b>	7.6 ng/mL
<b>Calibration Levels</b>	5-Point Calibration
<b>Reagent On-Board Stability</b>	<p>Unopened: ≥ 8 months when stored at 2-8°C</p> <p>Opened: 4 weeks when stored at 2-8°C</p>

## 25-OH VITAMIN D ASSAY PROCEDURE FOR CLINICAL CHEMISTRY ANALYZERS\*



\*Analyzer Dependent

Parameter questions for 25-OH Vitamin D Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email [support@diazyme.com](mailto:support@diazyme.com)

## ASSAY PRECISION

Precision was evaluated according to the CLSI EP5-A guideline. Controls and samples were measured daily over the span of 20 days, using three lots of reagents and one chemistry analyzer. 40 independent runs were performed on each specimen. Each run produced two measurements. 80 data points were obtained per specimen. Results are shown below:

25-OH Vitamin D (ng/mL)	Speciman	n	Within-run			Between-run		Total	
			Mean	SD	CV%	SD	CV%	SD	CV%
Control #1	80	23.1	1.47	6.4	1.04	4.5	1.68	7.3	
Control #2	80	45.7	2.06	4.5	1.67	3.7	2.12	4.6	
Sample #1	80	22.6	1.19	5.3	1.11	4.9	1.45	6.4	
Sample #2	80	31.7	1.42	4.5	1.59	5.0	1.81	5.7	
Sample #3	80	40.6	1.42	3.5	1.59	3.9	1.66	4.1	
Sample #4	80	48.6	2.32	4.8	1.71	3.5	2.41	4.9	
Sample #5	80	55.8	2.14	3.8	1.73	3.1	2.34	4.2	
Sample #6	80	65.4	2.03	3.1	1.79	2.7	2.42	3.7	
Sample #7	80	69.7	2.02	2.9	1.99	2.9	2.55	3.7	
Sample #8	80	92.8	2.52	2.7	2.02	2.2	3.40	3.7	
Sample #9	80	134.6	2.97	2.2	2.69	2.0	3.87	2.9	
Low Sample #1	80	9.4	1.22	13.0	0.98	10.4	1.31	14.0	
Low Sample #2	80	11.2	1.58	14.2	0.88	7.9	1.55	13.9	

## ASSAY INTERFERENCE

Interference studies were conducted according to the CLSI EP7-A2 guideline. The acceptance criterion was set at 10% or less deviation between the spiked sample and the control. The assay's results were not significantly affected by the following substances:

Substance	Concentration	Substance	Concentration
Conjugated Bilirubin	40 mg/dL	Ampicillin	1000 ng/mL
Free Bilirubin	40 mg/dL	Cyclosporine C	105 ng/mL
Hemoglobin	100 mg/dL	Cefoxitin	660 ng/mL
Ascorbic Acid	176 mg/dL	Acetylsalicylic Acid	1000 ng/mL
Triglycerides	750 mg/dL	Rifampicin	64 ng/mL
Uric Acid	20 mg/dL	Acetaminophen	200 ng/mL
Biotin	2 mg/dL	Ibuprofen	500 ng/mL
Human Serum Albumin	9 g/dL	Theophylline	100 ng/mL
N-Acetyl Cysteine Amide	1663 ng/mL		

Cross-reactivity of the Diazyme 25-OH Vitamin D Assay was determined by adding Vitamin D metabolites to serum pool samples. Based on the results in the table below, the assay did not cross react with Vitamin D2 and Vitamin D3 and the assay recovers both 25-OH Vitamin D2 and 25-OH Vitamin D3 similarly. Cross-reactivity with various Vitamin D metabolites is summarized in the table below:

Compound	Concentration tested	Cross-reactivity
25-OH Vitamin D3	44.0 ng/mL	100%
25-OH Vitamin D2	44.0 ng/mL	92.3%
Vitamin D3	44.0 ng/mL	1.0%
Vitamin D2	44.0 ng/mL	2.9%
1,25-(OH) <sub>2</sub> Vitamin D3	2.9 ng/mL	2.5%
1,25-(OH) <sub>2</sub> Vitamin D2	2.9 ng/mL	-1.5%
24R,25-(OH) <sub>2</sub> Vitamin D3	41.0 ng/mL	5.1%
3-epi-25-OH Vitamin D3	42.0 ng/mL	61.7%
3-epi-25-OH Vitamin D2	42.0 ng/mL	55.1%

\*%Cross-reactivity = (Corrected Assay Value/Concentration Spiked) \* 100

No significant cross-reactivity (4.1%) was found for Paricalcitol (Zemlar®) up to 25 ng/mL.

## DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064

PO Box 85608, San Diego, CA 92186

Tel: 858-455-4768 888-DIAZYME

[www.diazyme.com](http://www.diazyme.com) [sales@diazyme.com](mailto:sales@diazyme.com)



IGZ Instruments AG  
Furtbachstrasse 17  
8107 Buchs ZH

Tel: +41 44 456 33 33  
igz.ch igz@igz.ch

## DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland

Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366

[sales@diazyme.de](mailto:sales@diazyme.de)

## SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park

Shanghai, 201203, People's Republic of China

Tel: 086-21-51320668 Fax: 086-21-51320663

[www.lanyuanbio.com](http://www.lanyuanbio.com) [service@lanyuanbio.com](mailto:service@lanyuanbio.com)

