

CD34+ Control



CD-Chex CD34[®] is a positive procedural control for enumerating CD34 positive cells by flow cytometry. The CD34 positive cells have characteristics of human hematopoietic progenitor cells, including low to intermediate side scatter properties, CD34 antigen expression and low expression of CD45. CD-Chex CD34 is available in three levels providing reference values similar to levels found in bone marrow, cord blood and mobilized patient samples.

Compatibility

- BD Biosciences and Beckman Coulter[®] flow cytometry systems
- ISHAGE Guidelines, as well as commercially available enumeration kits

Features and Benefits

- Positive procedural control assayed for flow cytometry systems
- Ready to use; no dilution required
- Only commercially available control offered in three levels

Stability (Store 2 °C - 10 °C)

Closed-vial stability	90 days
Open-vial stability	30 days

STATS[®]

Interlaboratory Quality Control program features online accessibility with enhanced real-time reporting and improved account management capability. Open to all customers at no charge, the Streck *STATS* program includes personalized, easy-to-read reports that detail performance, identify trends and facilitate real-time prompt corrective action. More information can be found at www.streck.com.

Ordering Information

Level	Description	Catalog #
Level 1	1 x 1.0 mL	213336
Level 1	2 x 1.0 mL	213337
Level 2	1 x 1.0 mL	213346
Level 2	2 x 1.0 mL	213347
Level 1, Level 2	2 x 1.0 mL	213338
Level 1, Level 2	4 x 1.0 mL	213339
Level 3	1 x 1.0 mL	213382
Level 3	2 x 1.0 mL	213383



IGZ Instruments AG
Furlbachstrasse 17
8107 Buchs ZH

Tel. +41 44 456 33 33
igz.ch igz@igz.ch